COVID-19 Antigen Rapid Test Kit (Swab) For Self-Testing Package Insert

For self testing and in vitro diagnostic use only.

The COVID-19 Antigen Rapid Test Kit (Swab) is a single-use test kit intended to detect the novel coronavirus SARS-CoV-2 that causes COVID-19. This test is intended for home use with self-collected nasal swab samples in individuals aged 18 and older. Sampling from anyone under the age of 18 should be performed under the guidance of an adult. People who are unable to carry out the test on their own should seek support. This test utilizes the lateral flow immunoassay technology for the detection of nucleocapsid protein antigen in individuals with known or suspected COVID-19. Positive results are indicative of the presence of SARS-CoV-2. Individuals who test positive should self-isolate and seek additional care from their healthcare provider. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not rule out SARS-CoV-2 infection. Individuals who test negative and continue to experience COVID-like symptoms should seek follow up care from their healthcare provider.

PRINCIPLE

The COVID-19 Antigen Rapid Test Kit (Swab) is an immunochromatographic membrane assay that uses highly sensitive antibodies to detect SARS-CoV-2 nucleocapsid protein from nasal swab specimens. SARS-CoV-2 specific antibodies are immobilized onto the test region of the membrane and combined with other reagents/pads to construct a test strip. During testing, the specimen reacts with anti-COVID-19 antibodies conjugated to colored particles and precoated onto the sample pad of the test. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the reagents in test line region. Therefore, if the specimen contains COVID-19 antigen, a colored line will appear in test line. If the specimen does not contain COVID-19 antigen, no colored line will appear in the test line regions, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred. The test is designed to detect nucleocapsid protein antigens in nasal swab specimens, which is different from the mutation sites have occurred in the spike protein, so it is theoretically able to detect variants including those in UK, India, South Africa and Brazil.

KIT COMPONENTS

| Component | 1 Test /Kit | 5 Tests /Kit |
|-----------------------------|-------------|--------------|
| COVID-19 Antigen Test | 1 | 5 |
| Extraction tube with buffer | 1 | 5 |
| Sterilized nasal swab | 1 | 5 |
| Package insert | 1 | 1 |

ADDITIONAL SPECIAL EQUIPMENT

Timer

POTENTIAL RISKS AND BENEFITS

- Potential risks include:
 - Possible discomfort during sample collection.
 - Possible incorrect test results (see Read Results section).
- Potential benefits include:
 - The results, along with other information, can help your healthcare provider make informed recommendations about your care. The results of this test may help limit the spread of COVID-19 to your family and others in your community.

WARNINGS AND PRECAUTIONS

- Do not use after expiration date. Do not use if kit is damaged or open. Do not reuse the test.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Discard the used testing materials in accordance with local regulations.
- The extraction buffer contains a soft seal solution; if the solution contacts the skin or eyes, flush with copious amounts of water. Do not swallow the buffer. When swallowing the buffer, rinse the mouth thoroughly with water and give plenty of water to dilute the substance. If any discomfort, seek medical attention immediately.
- Children and elder please use the test under the guidance.

STORAGE AND STABILITY

Store unused test devices upon expiration at 4°C-30°C. If stored at 4°C-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

TEST PROCEDURE


Open the kit box, Check the components before use.

Please read all instructions carefully before you begin.

- Get a fat area ready, like a table. Make sure it is clear, clean and dry.
- Take off hand jewelry.
- Wash your hands for 20 seconds. Use soap and warm water, or hand sanitizer. Dry your hands using clean disposable paper towels.
- For better protection and avoid cross-contamination, disposable gloves, masks and eye protection (not provided in the package) are recommended.

- [Specimen Preparation]
- Press the hole on the box, peel off the foil film on the extraction tube and insert into the hole.
 - Open nasal swab package at the sticky end and take the nasal swab out.
 - Gently insert the soft tip of the nasal swab into both nostril about 2.5cm (1 inch) for adults. *Note: For child, the maximum depth of insertion into the nostril may be less than 2.5 cm and should be carefully and appropriately adjusted by the person, who collects sample.*
 - Firmly brush against the inside of the nostril in a circular motion 5 times or more.
 - Rotate the nasal swab to the right nostril and repeat the previous action. Make sure an adequate sample is collected.
 - Insert the nasal swab into the tube which contains extraction buffer.
 - Rotate nasal swab at least 8-10 times while pressing nasal swab tip against the bottom and the side of the tube.
 - Remove the nasal swab while squeezing and turning the nasal swab against the sides of the tube to release as much liquid as possible.
 - Cover the tube with dropper tip tightly and insert the tube back into the box.

- [Test Procedure]
- Open the sealed pouch and take out the test cassette. For best results, the test should be performed in one hour.
 - Hold the tube vertically upside down over the sample well.
 - Add 3 drops specimen into the sample well by gently squeezing the sides of the tube, then start the timer.

- Peel off the foil film on the extraction tube and insert into the hole.
- Insert nasal swab to an appropriate depth, not 5 times in each nostril.
- Rotate the swab 8-10 times.
- Press the swab and release as much liquid as possible.
- Place the cap.
- Add 3 drops to sample well(S).

- 4.Wait for colored lines to appear.The test result can be read in 10-15 minutes. DO NOT read after 20 minutes. [After the testing]
- After you have done the test, put all parts of the kit in the waste bag. Discard the waste bag in accordance with local regulations.
 - If you are doing more than 1 test, clean the table with 75% alcohol or sanitizer. Wash your hands between each test.

| POSITIVE RESULT | NEGATIVE RESULT | INVALID RESULT |
|---|---|--|
| <p>A colored line appears in the control line region (C) and a colored line appears in test line region (T). NOTE: The intensity of the color in the test line region will vary depending on the concentration of COVID-19 antigen in the specimen. Therefore, any shade of color in the test line region should be considered positive.</p> | <p>A colored line appears in the control line region (C) and no line appears in test line region (T).</p> | <p>No line appears in control line region (C). Insufficient buffer volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the procedure with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.</p> |

BUILT-IN CONTROL

This test contains a built-in control feature, the C line on the test. The C line develops after adding sample solution. Otherwise, review the whole procedure and repeat test with a new device.

WHAT SHOULD I DO AFTER TEST

| | |
|--------------------------------|---|
| If the test result is positive | There is currently a suspicion of a COVID-19 infection <ul style="list-style-type: none"> Contact your doctor / general practitioner or the local health department immediately Comply with local guidelines for self-isolation To have a PCR confirmatory test performed |
| If the test result is negative | Continue to comply with all applicable rules regarding contact with others and protective measures <ul style="list-style-type: none"> An infection may also be present if the test is negative If it is suspected, repeat the test after 1 - 2 days, as the coronavirus cannot be accurately detected in all phases of an infection |
| If the test result is invalid | Possibly caused by incorrect test execution <ul style="list-style-type: none"> Repeat the test If the test results remain invalid, contact a doctor or a COVID-19 test center |

Note: Do not take any decision of medical relevance without first consulting your medical practitioner

PERFORMANCE CHARACTERISTICS

- Limit of Detection: The Limit of Detection (LoD) was determined by evaluating different dilutions of heat inactivated SARS-CoV-2 virus. LoD is defined as the virus concentration - 1.3 x 10⁷ TCID₅₀/mL at which a minimum of 19 replicates out of 20 generate a Reactive result.
- High-dose Hook Effect: No hook effect was observed at 3.33X10⁷ TCID₅₀/mL of inactivated SARS-CoV-2 virus culture.
- Clinical study: A comparison was conducted with RT-PCR, according the test data listed below table, the Relative Sensitivity is 96.4% (212/220), the Relative Specificity is 100%(200/200) and the overall agreement is 98.1% ((212+200)/420).

| Safecare Test | PCR result | | Total |
|---------------|------------|----------|-------|
| | Positive | Negative | |
| Positive | 212 | 0 | 212 |
| Negative | 8 | 200 | 208 |
| Total | 220 | 200 | 420 |

4. Cross-reactivity: Cross-reactivity studies were performed to demonstrate that the test does not react with the following microorganisms in the table below at concentration of 1x10⁷ TCID₅₀/mL for viruses and 1x10⁸ CFU/mL for bacteria.

| | | | | |
|------------------------------|-----------------------------|-------------|-------------------------------|-----------------------------------|
| Human metapneumovirus (hMPV) | Human parainfluenza virus 1 | Adenovirus | <i>Bordetella pertussis</i> | <i>Streptococcus pneumoniae</i> |
| Human coronavirus OC43 | Human parainfluenza virus 2 | Rhinovirus | <i>Chlamydia pneumoniae</i> | <i>Streptococcus pyogenes</i> |
| Human coronavirus 229E | Human parainfluenza virus 3 | Enterovirus | <i>Haemophilus influenzae</i> | <i>Mycobacterium tuberculosis</i> |
| Human coronavirus NL63 | Human parainfluenza virus 4 | Influenza A | <i>Legionella pneumophila</i> | <i>Staphylococcus aureus</i> |
| Respiratory Syncytial Virus | MERS | Influenza B | <i>Mycoplasma pneumoniae</i> | <i>Candida albicans</i> |

- Interference: The following endogenous interference substances were evaluated at the concentrations listed and no effect was found.

Whole blood (2%), three OTC nasal sprays (10%), three OTC nasal drops (25%), three nasal mouthwashes (25%), 4-Acetamidophenol (10mg/mL), Acetylsalicylic acid (20mg/mL), Chlorpheniramine (5 mg/mL), Dextromethorphan (10mg/mL), Diphenhydramine (5mg/mL), Ephedrine (20mg/mL), Guaiacol glyceryl ether (20mg/mL), Oxymetazoline (100mg/mL), Phenylephrine (100mg/mL), Pseudoephedrine (20mg/mL), Oseltamivir Phosphate (10mg/mL), Mupirocin (10mg/mL), Vitamin A (10%), D-Panthenol (10%)

LIMITATIONS AND POSSIBLE ERRORS

- The COVID-19 Antigen Rapid Test Kit (Swab) is intended for use as a self-test and may only be used for qualitative detection of SARS-CoV-2 antigens. The colour intensity of a positive line shall not be evaluated as quantitative or semi-quantitative.
- The COVID-19 Antigen Rapid Test Kit (Swab) should be only for the detection of SARS-CoV-2 antigen, not for any other viruses or pathogens.
- The performance was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.
- A negative result does not exclude the possibility of COVID-19 infection.
- The results obtained with this assay, especially in the case of weak test lines which are difficult to interpret, should retest or go to a medical institution for testing.
- The test is intended for infection detection and not for determine infection status.The test is used for the auxiliary diagnosis of patients with the COVID-19 and can not be used as the sole diagnostic indicator of whether the test subject is infected with the COVID-19.

ACCESSORIES

| Accessory | Manufacturer | EC Representative | CE Mark |
|-----------|---|---|---------|
| Swab A | Jiangsu Changfeng Medical Industry Co.,Ltd. Touqiao Town,Guangding District Yangzhou Jiangsu 225109, P.R.China | Lins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany | |
| Swab B | Goodwood Medical Care Ltd., 1-2 Floor,3-619 Yongzheng Street ,Jinzhou District,Dalian 116100 Liaoning, P.R.China | CMC Medical Devices & Drugs S.L., C/Choracio Lengo No.18,CP29006,Málaga,Spain | |
| Swab C | QITO TEST LABWARE MANUFACTURING CO.,Ltd. No.339 Beihai West Road,Haimen, 226100 Jiangsu, P.R. China | Wellking Ltd Enterprise Hub,NW Business Complex,1 Beraghmore Road,Derry,BT48 8SE Northern Ireland | |
| Swab D | Jiangsu Hanteng Medical Technology Co.,Ltd. 16-B4,#1 North Qingyang Road, Tiaaming District, Changzhou, 213017 Jiangsu P.R. China | Luxus lebenswelt GmbH Koshtstr.1,47877,Willich,Germany | |
| Swab E | Jiangsu Rongye Technology Co.,Ltd Touqiao Town,Yangzhou City,225109 Jiangsu, P.R.China | Riomavix S.L. Calle de Almansa 55,1D, Madrid 28039 Spain | |

INDEX OF SYMBOLS

| | | | |
|--|---|--|--|
| | Do not reuse | | In vitro diagnostic use only |
| | Stored between 4-30°C | | Consult instruction for use |
| | Keep away from sunlight | | Lot number |
| | Use-by date | | Contains sufficient for <-> tests |
| | Manufacturer | | Date of Manufacture |
| | Sterilized using ethylene oxide | | Sterilized using irradiation |
| | Authorized Representative in the European Community | | CE Mark |
| | Importer | | Safecare Biotech(Hangzhou) Co., Ltd. Building 2/203, No.18 Haisu Rd., Guangning Sub-district Yuhang District, Hangzhou, 311121, China Tel/Fax: +86 571 81389219 Email: admin@safecare.com.cn |
| | NIC GmbH | | EC REP |
| | Angela Q. GmbH | | EC REP |

| | | | | |
|--|--|--|----------|---|
| | Safecare Biotech(Hangzhou) Co., Ltd. Building 2/203, No.18 Haisu Rd., Guangning Sub-district Yuhang District, Hangzhou, 311121, China Tel/Fax: +86 571 81389219 Email: admin@safecare.com.cn | | NIC GmbH | Erkerweg 13, 49076 Osnabrück, Germany Tel:+49 541 9116706 Email:info@nic-industry.com |
| | Angela Q. GmbH | | REP | EC REP |
| | Angela Q. GmbH | | REP | EC REP |

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



ACCESSOIRES

| Accessoire | Fabricant | Représentant autorisé dans la Communauté européenne | Marquage CE |
|--------------|---|--|---|
| Écouvillon A | Jiangsu Changfeng Medical Industry Co.,Ltd. Toujiao Town,Guangling District Yangzhou Jiangsu 225109, P.R.China | lins Service & Consulting GmbH Obere Seegasse 34/2,69124 Heidelberg Germany |  |
| Écouvillon B | Goodwood Medical Care Ltd, 1-2 Floor,3-919 Yongzheng Street,Jinzhou District,Dalian 116100 Liaoning, P.R.China | CMC Medical Devices & Drugs S.L. CHirolaco Lengo No.18,CP29006,Málaga, Spain |  |
| Écouvillon C | CITOTEST LABWARE MANUFACTURING CO.,Ltd. No.339 Beihai West Road,Haimen, 226100 Jiangsu, P.R. China | WellKang Ltd Enterprise Hub,NW Business Complex.1 Berghmore Road,Derry,BT48 8SE Northern Ireland |  |
| Écouvillon D | Jiangsu Han-Heng Medical Technology Co.Ltd, 16-B4-#1 North Qingyang Road, Tianjing District, Changzhou, 213017 Jiangsu P.R. China | Luxus lebenswelt GmbH Koshstr.1,47877,Willich,Germany |  |
| Écouvillon E | Jiangsu Rongye Technology Co.,Ltd Toujiao Town,Yangzhou City,225109 Jiangsu, P.R.China | Riomaviv S.L. Calle de Almansa 55,1D, Madrid 28039 Spain |  |

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